

ClassicStar® NIV Mask Series: Nasal masks, vented and non-vented

More comfort with anatomical fit and effective seal.



The ClassicStar NIV nasal masks are new additions to the existing portfolio of ClassicStar NIV full-face masks.

Full-face masks are often the first-line strategy in the initial management of acute respiratory failure as they cover both the nose and the mouth, so the concerns relating to mouth breathing are kept to a minimum. However, nasal masks are often chosen if non-invasive ventilation has to be prolonged due to the improved comfort and reduced complications associated with full-face masks. Nasal masks are mostly preferred for patients where mouth breathing is not a significant problem.¹⁾ They may also provide valuable support for non-invasive ventilation therapy in the home care environment.


The ClassicStar NIV nasal masks are for single-patient use. They are available in 3 sizes.

The non-vented mask is designed for ventilators which have an integral safety valve (active exhaust system). The non-vented mask is suitable for the following Dräger ventilators: Evita Series, Savina (both with NIV Option) and Carina. The vented mask includes an exhalation port to provide adequate flow to wash out exhaled CO₂

and is suitable for ventilators with continuous positive airway pressure without active exhalation system.

- The ClassicStar NIV nasal mask has an air-filled mask cushion which can be adjusted individually to the patient's face. The mask cushion can be inflated or deflated by means of the Dräger pump ball. Due to its adaptability, the ClassicStar NIV nasal mask provides increased comfort and effective seal (see p. 2).
- The forehead support and forehead pad are adjustable to match the patient's forehead curve and thus help maximizing patient comfort and mask seal.
- The 4-point headgear provides a stable and secure fixation reducing mask movement or displacement.
- The ClassicStar mask swivel is moveable through 360°. In addition, the swivel can be moved vertically providing more flexibility.

¹⁾ Girault, Christophe, et al.: Interface Strategy during noninvasive positive pressure ventilation for hypercapnic acute respiratory failure, in: Crit. Care Med., 2009, 37(1): 124-131



ClassicStar® NIV nasal mask non-vented

ClassicStar® NIV nasal mask vented



ClassicStar® NIV Full-Face Mask with Standard Elbow (SE).



Two-sided pump ball for air-management.

The Dräger Sizing Gauge helps you to choose the right mask size – there is no need to open the mask package.

Finding the right mask size takes only 2 steps: First, the length of the patient's nose should be measured by holding the sizing gauge against the side of the patient's nose. Second, by simply turning the sizing gauge clockwise and holding it against the bottom of the patient's nose,

you can measure the width of the patient's nose.

The colors follow the well-known traffic light principle: red for small, yellow for medium, and green for large sizes. Depending on which color the lengths and width of the patient's nose match, you can pick a small, medium or large nasal mask. The masks are packaged individually, fully assembled, and ready to be used.



D-9952-2009

Individual adaptation (formability) of ClassicStar® NIV nasal masks, vented and non-vented



D-12187-2009



D-12191-2009



D-12182-2009

The mask can be adjusted to the patient's face by increasing or decreasing the amount of air in the mask cushion. After adjusting, the mask stays in position.

ORDER INFORMATION

Description	Packing unit	Order number
ClassicStar non-vented nasal mask, for non-invasive ventilation, size S	1 pc.	MP01623
ClassicStar non-vented nasal mask, for non-invasive ventilation, size M	1 pc.	MP01624
ClassicStar non-vented nasal mask, for non-invasive ventilation, size L	1 pc.	MP01625
ClassicStar vented nasal mask, for non-invasive ventilation, size S	1 pc.	MP01620
ClassicStar vented nasal mask, for non-invasive ventilation, size M	1 pc.	MP01621
ClassicStar vented nasal mask, for non-invasive ventilation, size L	1 pc.	MP01622
Two-sided pump ball for air-management	1 pc.	MP01590
Sizing Gauge for ClassicStar NIV nasal masks vented and non-vented	3 pcs.	MP01632

HEADQUARTERS
 Drägerwerk AG & Co. KGaA
 Moislinger Allee 53–55
 23558 Lübeck, Germany

 www.draeger.com

REGION EUROPE CENTRAL AND EUROPE NORTH
 Dräger Medical GmbH
 Moislinger Allee 53–55
 23558 Lübeck, Germany
 Tel +49 451 882 0
 Fax +49 451 882 2080
 info@draeger.com

REGION MIDDLE EAST, AFRICA, CENTRAL AND SOUTH AMERICA
 Dräger Medical GmbH
 Branch Office Dubai
 Dubai Healthcare City, P.O. Box 505108
 Dubai, United Arab Emirates
 Tel + 971 436 24 762
 Fax + 971 436 24 761
 contactuae@draeger.com

REGION NORTH AMERICA
 Draeger Medical, Inc.
 3135 Quarry Road
 Telford, PA 18969-1042, USA
 Tel +1 215 721 5400
 Toll-free +1 800 437 2437
 Fax +1 215 723 5935
 info.usa@draeger.com

REGION EUROPE SOUTH
 Dräger Médical S.A.S.
 Parc de Haute
 Technologie d'Antony 2
 25, rue Georges Besse
 92182 Antony Cedex, France
 Tel +33 1 46 11 56 00
 Fax +33 1 40 96 97 20
 dlmfr-contact@draeger.com

REGION ASIA / PACIFIC
 Draeger Medical
 South East Asia Pte Ltd
 25 International Business Park
 #04-27/29 German Centre
 Singapore 609916, Singapore
 Tel +65 6572 4388
 Fax +65 6572 4399
 asia.pacific@draeger.com

Manufacturer:
 Dräger Medical GmbH
 23542 Lübeck, Germany
 The quality management system at Dräger Medical GmbH is certified according to ISO 13485, ISO 9001 and Annex II.3 of Directive 93/42/EEC (Medical devices).