

BabyFlow™ nasal CPAP for Babylog 8000

Easy application and
stable CPAP

The BabyFlow® nasal CPAP system for the Babylog 8000 ventilator series allows easy application and stable CPAP at low operating costs due to the use of standard patient tubing and reusable parts. The unique Prong design and attachment accessories are designed to reduce trauma associated with the delivery of nasal CPAP

BabyFlow Features:

- Unique nasal prong design
- Unique headbands
- Innovative attachment technique





BabyFlow™ Features

- Seamless transition to nCPAP with Babylog 8000 series ventilators
- Low noise level during nCPAP application supports developmental care
- Flow independent nCPAP control
- Minimal deadspace
- No need for separate nCPAP device
- Helps to reduce equipment complexity at the bedside
- Inventory is simplified and costs reduced
- Easy training and operation

Unique headbands

- A complete line of six color coded headbands available
- Headband allows easy access during head ultrasound and IV cannulation procedures
- Each headband has 7 separate velcro attachment points and minimal seams for max. flexibility and patient comfort
- Headbands, all handmade from soft microfiber material (91 % microfiber, 9 % Elastan)
 - Microfiber is resistant to humidity
- All headbands are washable at 60 °C

Innovative attachment technique

- Attachment design adjusts the BabyFlow with minimal disturbance to the baby
- Standard delivery circuit mounts on the headband, allows for free movement of the baby and for Kangaroo Care
- Seven velcro attachment points are built into the headband to allow for quick and easy securing of the BabyFlow to the infant
- Prongs remain in place more effectively than with other systems



Unique anatomical prong and mask design

- Nasal silicone prongs in three different sizes (S,M,L) allows for a perfect fit
- Greater comfort for your patient due to optimal anatomical design

Headband sizing chart

Size	Color	Head circumference	Approx Weight in grams
XS	White	21 cm	500 - 700
S	Yellow	23 cm	700 - 1000
M	Red	25 cm	1000 - 1300
L	Green	28 cm	1300 - 1600
XL	Orange	31 cm	1600 - 2400
XXL	Blue	34 cm	2400 - 3000

Order information

Name and description	Order no.
Ask for BabyFlow Starterpackage	84 18 492
Consisting of:	
2 x BabyFlow circuits	84 18 470
2 x Prongs small (packs of 10 pcs.)	84 18 415
2 x Prongs medium (packs of 10 pcs.)	84 18 416
2 x Prongs large (packs of 10 pcs.)	84 18 417
2 x Neo masks (packs of 10 pcs.) > 700 - ca. 1500 gr.	84 18 491
1 x Headband XS	84 18 418
1 x Headband S	84 18 419
3 x Headband M	84 18 420
2 x Headband L	84 18 421
1 x Headband XL	84 18 422
1 x Headband XXL	84 18 423
Each component of the starter package is also available separately	

Europe, Middle East, Africa:

Dräger Medical AG & Co. KGaA
 Moislinger Allee 53-55
 23542 Lübeck
 GERMANY
 Phone: +49-1805-3 72 34 37
 Fax: +49-451-882-37 79
 E-mail:
 Business.Support@draeger.com

North America:

Draeger Medical, Inc.
 3135 Quarry Road
 Telford, PA 18969
 USA
 Phone: +1-215-721-5400
 Toll-free: +1-800-437-2437
 Fax: +1-215-723-5935
 E-mail:
 wwwinfo@draegermed.com

Spain / Latin America:

Dräger Medical Hispania S.A.
 c/ Xaudaró n° 5
 28034 Madrid
 SPAIN
 Phone: +34-91-728 34 00
 Fax: +34-91-358 51 61
 E-mail Spain:
 draegermedical@draeger.es
 E-mail Latin America:
 latinoamerica@draeger.es

Asia / Pacific:

Draeger Medical Asia Pacific Ltd.
 Unit 2205, 22nd floor
 Harbour Centre
 25 Harbour Road
 Wanchai
 HONG KONG
 Phone: +852-28 77 30 28
 Fax: +852-28 77 33 37
 E-mail: apinfo@draeger.com.hk