

Quickstop hinged arms

Keeping everything
close at hand

To leave your hands free to care for the patient, all the joints of the Quickstop hinged arm can be fixed into place and released with just one quick hand movement. The hinged arm is totally flexible, allowing it to be positioned precisely where it is needed. All cables, lines and hose systems can be arranged in such a way as to more or less completely prevent accidental disconnection, making the Quickstop hinged arm an essential accessory in anaesthesia, intensive care, surgery or indeed just about any other area of the hospital. The arm is provided with a universal rail bracket which can be used to attach it to all kinds of different objects (operating tables, beds, wall-mounted rails and devices). Thanks to the high quality of its materials, you will appreciate the benefits of this durable and maintenance-free Dräger product for years to come. Quickstop hinged arms - keeping everything close at hand.

- Quick
- Precise
- Flexible
- Universal
- Robust



Quickstop hinged arm 1

Length 1067 mm, contains:
3 x hose clamps (single)
3 x hose clamps (double)
Order no. 2M 85 705



Quickstop hinged arm 2

Length 1078 mm, contains:
2 x hose clamps (double)
1 x hose clamp with ball
Order no. 2M 85 706

Optional accessories
Holder for Y-piece
Order no. 84 09 746



Quickstop hinged arm, pressure sensor

Length 463 mm
Order no. 2M 86 464



Quickstop hinged arm, short

Length 403 mm, contains:
1 x universal tube holder plate for
several cable- and tube dimensions
Order no. 2M 86 463



Spare Parts

Hose clamp (single)
Order no. 84 08 356
Hose clamp (double)
Order no. 84 06 138
Hose clamp with ball
Order no. 84 09 841

Europe, Middle East, Africa, Latin America,
Asia, Pacific:

Dräger Medical AG & Co. KGaA
Moislinger Allee 53-55
23542 Lübeck
GERMANY
Tel: +49-1805-3 72 34 37
+49-451-882-808
Fax: +49-451-882-37 79
E-mail:
Business.Support@draeger.com

Manufacturer: Dräger Medical AG & Co. KGaA,
23542 Lübeck, Germany

The quality management system at
Dräger Medical AG & Co. KGaA is certified
according to ISO 13485, ISO 9001 and
Annex II of Directive 93/42/EEC (Medical devices).

www.draeger-medical.com